# FOOT ORTHOSES

# INSTRUCTIONS FOR USE

END CONSUMER / PATIENT

ES – PLANTILLAS ORTOPÉDICAS – INSTRUCCIONES DE USO – CLIENTE FINAL  
NL – ORTHOPEDISCHE INLEGZOLEN – GEBRUIKSAANWIJZING – EINDGEBRUIKER  
FR – ORTHÈSES PLANTAIRES – INSTRUCTIONS D'UTILISATION – CLIENT FINAL  
SE – INLÄGG – BRUKSANVISNING – SLUTKUND  
DE – EINLAGEN – GEBRAUCHSANWEISUNG – ENDKUNDE

DEUTSCH

SVENSKA

FRANÇAIS

NEDERLANDS

ESPAÑOL

ENGLISH

S.54

S.44

S.34

S.24

S.14

S.4

## CONTENTS

Introduction .....	6
Product description .....	7
Purpose .....	7
Indications .....	7
Contraindications .....	8
Instructions for use and information on risks .....	8
Safety instructions and residual risks .....	10
Cleaning and care .....	11
General rules .....	11
Technical specifications / Parameters .....	12
Disposal .....	12
Guarantee .....	13

**For visually impaired people, this document can be downloaded in a pdf.-format on [www.perpedes.com](http://www.perpedes.com).**

## INTRODUCTION

Dear clients,

Please read the instructions for use carefully before using the product. Follow all instructions, especially the safety instructions.

Perpedes GmbH  
Haertwasen 8-14  
73252 Lenningen  
Germany

T. +49 (0)7021 738 30-0  
info@perpedes.com  
www.perpedes.com

## PRODUCT DESCRIPTION

Foot orthotics are to be used on the sole of the feet and are best be worn with footwear that is indicated as suitable for use with orthotics or „insole-friendly“. The adjustment to the foot and the correct fit of the orthotic's functional elements must be carried out by trained specialist personnel. The orthotics must also be fitted in the shoes. The outline of an orthosis corresponds to the shape of the shoe in which it is are worn in.

## PURPOSE

After user-specific fitting and customization, foot orthotics serve to prevent or aid in the conservative treatment of foot deformities or complaints. They are worn in footwear suitable for use with orthotics. They exercise a supporting, guiding or muscle-stimulating effect to preserve the arches of the foot. Foot orthoses also relieve specific sections of the foot and alleviate complaints from overuse.

## INDICATIONS

Indications for the use of foot orthoses include deformities and complaints that affect feet and lower limbs and require support of these structures or relief of certain areas of the feet.

## CONTRAINDICATIONS

Contraindications to the use of foot orthoses include:

- Deformities of the feet with strong variations in foot proportion and shape
- Distinct sensory disturbance of the feet

## INSTRUCTIONS FOR USE AND INFORMATION ON RISKS

Before being used for the first time, the foot orthotics must be fitted to the needs of the user/patient by qualified personnel.

- Foot orthoses are to be worn with (wrinkle-free) socks, on the correct side (left orthotic on the left foot, right orthotic on the right foot) and always in pairs! The orthotics must not slip or buckle when putting on the shoes.
- The shoe must provide sufficient internal space (length/width, toe space [foot length + 1 cm ]) for both the foot and the orthotic. Any inlays or soft layers that come with the shoe must be removed before inserting foot orthotics to ensure a correct, plane contact surface and sufficient space available.
- The forefoot, midfoot and hindfoot must be supported by the orthotic across the entire surface and must not hang over the edges; otherwise, blistering or callus formation may occur on the soles or sides of the feet.
- Materials used for foot orthoses are subject aging or wear and tear and may change in their characteristics. Foot

orthoses must not be worn if the covers, padding material or orthotic shell have become perforated or wrinkled or broken due to overuse!

- Foot orthotics are influenced by body weight, pressure, torsion and friction and therefore have a load-dependent lifespan. Temperatures, hosiery materials, cover materials, footwear and environmental influences also limit their service life.
- It is recommended that foot orthoses are to be checked after 6 months at the latest for function and condition by an orthopedic technician, podiatrist, or other qualified specialized personnel.
- Feet should be checked regularly for signs of pressure, redness and other changes.

## SAFETY INSTRUCTIONS AND RESIDUAL RISKS

Fitting, trimming and modifications of foot orthoses as instructed in this manual should only be carried out by qualified personnel.

The product is only for use with one patient/user.

It should be ensured that the foot orthoses are worn side-correct and in footwear with sufficient internal space. Foot orthoses may not be worn if they are wrinkled or have slipped out of place in the shoe.

A break-in period is normal, but if excessive pain, swelling, numbness or changes in sensitivity, unusual reactions or discomfort occur during use, please contact your physician, orthopedic technician/podiatrist or qualified personnel immediately.

The product must not be used in water and must be protected from heat  $> 50^{\circ}\text{C}/122^{\circ}\text{F}$ .

If the product is damaged, please replace it immediately or have it repaired by qualified personnel.

Depending on the country, foot orthoses may be medically prescribed products that should be worn under supervision by a physician or other qualified personnel. However, orthotics may only be used in accordance with these instructions for use and within their intended purpose.

## CLEANING AND CARE

- Foot orthotics need to be taken out of the shoes daily to allow moisture to evaporate.
- If necessary, the top cover materials can be cleaned with mild washing lotion and a damp sponge or cloth. The soap residue should be removed. Do not use external heat sources to dry; only allow to air-dry. Orthotics are not suitable for washing machines!
- Mild disinfectants (suitable for human skin) can be sprayed on/applied sparingly. Discoloration of the cover material may occur and does not constitute a defect under warranty.

## GENERAL RULES

- Foot orthoses may only be used by one patient/user due to indication-based selection, individual fitting and hygienic reasons. Orthotics are therefore not suitable for multiple use (re-use for different people).
- Check the product regularly, but at least every 4 weeks for wear and tear. The properties of orthotics diminish as a result of repeated use (wearing). It is therefore not possible to define a specific lifespan or wearing period. Foot orthoses should be replaced when worn.
- The combination of foot orthoses with other medical devices is only allowed with selected and adaptable products. Maintenance must be carried out by qualified specialist personnel!

## TECHNICAL SPECIFICATIONS / PARAMETERS

Foot orthotic materials may vary depending on the model. However, a foot orthosis is composed of the following main components:

- Self-supporting orthotic shell or orthotic base made of EVA\* or cork
- Padding materials made of EVA\*, PUR\* or SBR\* (optional)
- Reinforcement materials (optional)
- Supporting components, e.g. met pad made of foam rubber (optional)
- Top and base covers made of PUR\*, PA\* or PES\*
- Posture-affecting "Balancer" heel module made of TPU\* (optional)

\*EVA Ethylene-Vinyl Acetate

\*PUR Polyurethane

\*SBR Styrene Butadiene Rubber

\*PA Polyamide

\*PES Polyester

\*TPU Thermoplastic Polyurethane

## DISPOSAL

You may dispose of the orthotics in accordance with the applicable national regulations.

## GUARANTEE

The current statutory provisions apply. Disregard of the instructions for use, improperly performed cleaning as well as independent technical changes and additions to the product will void the warranty and the general product liability.

The CE marking confirms compliance with the requirements of Regulation (EU) 2017/745 on medical devices (MDR).

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G.L.P. Graupis Leather Products Ltd.
c/o Perpedes Canada Inc.